



## Rabbit anti-T-bet

**Cat. No.: AIB-30191 (0.5 ml Concentrate);  
AIB-30235 (6 ml Ready-to-use)**

### Instructions for use

#### Intended use

This antibody is designed for the specific localisation of T-bet protein in formalin-fixed, paraffin-embedded tissue sections and in frozen sections.

Anti-T-bet antibody is intended for in vitro diagnostic use.

#### Specifications

<b>Specificity:</b>	T-bet
<b>Clone:</b>	ZSTB
<b>Isotype:</b>	Rabbit IgG1
<b>Species reactivity:</b>	Human +, others not tested

#### Summary and Description

T-bet is a T-box transcription factor and a marker for mature T-cells. The protein is expressed in very small amounts in Thp (T-helper precursor) cells and is absent in T-lymphoblastic leukaemia and lymphoma cells. In lymphocytes of the interfollicular T-cell zone reactive lymphoid tissues (tonsil, lymph nodes, and spleen) show a nuclear staining. Germinal centres, mantle and marginal zones are negative for T-bet.

T-bet can be found in some B-cell lymphoproliferative diseases, especially in early phases of B-cell development. It is also found in neoplasms derived from mature B-cells like CLL/SLL, marginal zone lymphoma and hairy cell leukaemia. On the other hand B-cell neoplasms derived from pre-germinal centre B-cells or germinal centre B-cells like mantle cell lymphoma, follicular lymphoma, diffuse large B-cell lymphoma and Burkitt lymphoma are negative for T-bet. Therefore antibodies against T-bet can be useful tools for diagnosis and sub-typing of B- and T-cell proliferative diseases.

#### Reagent provided








Rabbit monoclonal antibody from tissue culture supernatant in PBS, pH 7.4 including carrier protein and preservative for stabilisation in the following formats:

<b>Concentrate:</b>	0.5 ml	(Cat. No. AIB-30191)
<b>Ready-to-use:</b>	6 ml	(Cat. No. AIB-30235)

#### Dilution of primary antibody

Dilution of Nordic Biosite' concentrated antibody depends on the detection system used. The final working dilution must always be determined by the user. The elaboration of staining protocol should be done by an experienced specialist. For Nordic Biosite' recommendations see chapter 'Staining procedure'.

Explanations of the symbols on the product label:

	Catalog Number Bestellnummer Reference du catalogue		Batch Code Chargenbezeichnung Code du lot	 <b>Manufacturer</b>  Nordic BioSite AB Propellervägen 4A S-183 62 Täby Sweden Tel: +46 (0)8 5444 33 40 Fax: +46 (0)8 756 94 90 info@nordicbiosite.com www.nordicbiosite.com
	Use By Verwendbar bis Utiliser jusqu'à		In Vitro Diagnostic Medical Device In vitro Diagnostikum Dispositif médical de diagnostic in vitro	
	Consult Instructions for use Gebrauchsanweisung beachten Consulter les instructions d'utilisation		Temperature Limitation Lagerungstemperatur Limites de température	

## Storage and handling

The antibody should be stored at 2-8°C without further dilution.

Dilutions of the concentrated antibody should be done in a suitable antibody dilution buffer (e.g. BCB-20005/BCB-20006 from Nordic Biosite). The diluted antibody should be stored at 2-8°C after use. The stability of this working solution depends on various parameters and has to be confirmed by appropriate controls.

The antibody provided is suitable for use until the expiry date indicated on the label, if stored at 2-8°C. Do not use product after the expiry date. Positive and negative controls should be run simultaneously with all specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Nordic Biosite's technical support or your local distributor.

## Precautions

Use through qualified personnel only.

Wear protective clothing to avoid contact of reagents and specimens with eye, skin and mucous membranes. If reagents or specimens come in contact with sensitive area, wash with large amounts of water.

Microbial contamination of the reagent must be avoided, since otherwise non-specific staining may occur.

Sodium azide (NaN<sub>3</sub>), used for stabilisation, is not considered hazardous material in the concentration used.

Reaction of sodium azide with lead or copper in drainage pipes can result in the formation of highly explosive metallic azides. Sodium azide should be discarded in a large volume of running water to avoid formation of deposits. A material safety data sheet (MSDS) for the pure substance is available upon request.

## Staining procedure

Refer to the following table for conditions specifically recommended for this antibody. Also refer to detection system data sheets for guidance on specific staining protocols or other requirements.

### Parameters

\*Pre-treatment

20019/-20020)

\*Control tissue

\*Working dilution

\*Incubation time

### Nordic BioSites recommendations

Heat Induced Epitope Retrieval (for example in T-EDTA Buffer pH 9.0 (BCB-

Tonsil, hairy cell leukaemia

1:100-1:500 (for concentrates)

30-60 minutes

## Quality control

The recommended positive control tissues for this antibody are tonsil or hairy cell leukaemia. We recommend carrying out a positive and a negative control with every staining run. Please refer to the instructions of the detection system for guidance on general quality control procedures.

## Troubleshooting

If you observe unusual staining or other deviations from the expected results please read these instructions carefully, refer to the instructions of the detection system for relevant information or contact your local distributor.








## Expected results

This antibody stains positive in the nuclei of positive cells in formalin-fixed, paraffin-embedded tissue sections.

Further details about the expression pattern of T-bet can be found in the chapter 'Summary and Description'.

Interpretation of the staining results is solely the responsibility of the user. Any experimental result should be confirmed by a medically established diagnostic procedure.

Explanations of the symbols on the product label:

	Catalog Number Bestellnummer Reference du catalogue		Batch Code Chargenbezeichnung Code du lot	 Manufacturer  Nordic BioSite AB Propellervägen 4A S-183 62 Täby Sweden Tel: +46 (0)8 5444 33 40 Fax: +46 (0)8 756 94 90 info@nordicbiosite.com www.nordicbiosite.com
	Use By Verwendbar bis Utiliser jusqu'à		In Vitro Diagnostic Medical Device In vitro Diagnostikum Dispositif médical de diagnostic in vitro	
	Consult Instructions for use Gebrauchsanweisung beachten Consulter les instructions d'utilisation		Temperature Limitation Lagerungstemperatur Limites de température	

## Limitations of the Procedure

Immunohistochemistry is a complex technique involving both histological and immunological detection methods. Tissue processing and handling prior to immunostaining, for example variations in fixation and embedding or the inherent nature of the tissue can cause inconsistent results (Nadji and Morales, 1983). Endogenous peroxidase, alkaline phosphatase or biotin may cause non-specific staining depending on the detection system used. Tissues containing Hepatitis B Surface Antigen (HBsAg) may give false positive results with HRP (horse radish peroxidase) detection systems (Omata et al, 1980). Inadequate counterstaining and mounting can influence the interpretation of the results.

Nordic Biosite warrants that the product will meet all requirements described from its shipping date until the expiry date is reached, if the product is stored and utilised as recommended. No additional guarantees can be given. Under no circumstances shall Nordic Biosite be liable for any damages arising out of the use of the reagent provided.

## Performance characteristics

Nordic Biosite has conducted studies to evaluate the performance of the antibody for use with a standard detection system. The product has been found to be sensitive and specific to the antigen of interest with minimal or no cross-reactivity.

## Bibliography

Jöhrens K et al. Am J Surg Pathol 31:1181-5, 2007  
Atayar C et al. Am J Pathol 166:127-34, 2005  
Harashima A et al. Leuk Res 29:841-8, 2005  
Dorfman DM et al. Am J Clin Pathol 122:292-7, 2004  
Marafioti T et al. Am J Pathol 162:861-71, 2003








Zhang WX, Yang SY. Genomics 70:41-8, 2000  
Szabo SJ et al. Cell 100:665-69, 2000  
Nadji M and Morales AR Ann N.Y. Acad Sci 420:134-9, 1983  
Omata M et al. Am J Clin Pathol 73(5): 626-32, 1980

November 11, 2013

Rev: A1113

Doc: DBE\_AIB-30191/-30235

Explanations of the symbols on the product label:

	Catalog Number Bestellnummer Reference du catalogue		Batch Code Chargenbezeichnung Code du lot	 Manufacturer  Nordic BioSite AB Propellervägen 4A S-183 62 Täby Sweden Tel: +46 (0)8 5444 33 40 Fax: +46 (0)8 756 94 90 info@nordicbiosite.com www.nordicbiosite.com
	Use By Verwendbar bis Utiliser jusqu'à		In Vitro Diagnostic Medical Device In vitro Diagnostikum Dispositif médical de diagnostic in vitro	
	Consult Instructions for use Gebrauchsanweisung beachten Consulter les instructions d'utilisation		Temperature Limitation Lagerungstemperatur Limites de température	