



## Mouse anti-CD138 (Syndecan-1)

**Cat. No.: AIB-30018 (1 ml Concentrate); AIB-30019 (0.5 ml Concentrate);  
AIB-30017 (6 ml Ready-to-use)**

### Instructions for use

#### Intended use

This antibody is designed for the specific localisation of Syndecan-1 (CD138) in formalin-fixed, paraffin-embedded tissue sections.

Anti-CD138 antibody is intended for in vitro diagnostic use.

#### Specifications

<b>Specificity:</b>	CD138 (Syndecan-1)
<b>Clone:</b>	B-A38
<b>Isotype:</b>	Mouse IgG1
<b>Species reactivity:</b>	Human +, others not tested

#### Summary and Description

CD138 (Syndecan-1) is a trans-membranous proteoglycan which is detectable on the cytoplasmic membrane of nearly 95% of plasma cells in formalin-fixed paraffin-embedded tissue sections. Other haematopoietic cells, endothelial cells and lymphoplasmacytoid lymphomas are negative for CD138. CD138 is expressed by pre-B-cells. It is not detectable during further B-cell differentiation and reappears again in the plasma cell stadium (Sanderson *et al.* 1989).

This anti-CD138 antibody is a useful tool for labelling of plasma cells and characterisation of plasmacytoid lymphomas. Various forms of Hodgkin's disease have also been shown positive staining with this antibody.

#### Reagent provided

Mouse monoclonal antibody in PBS with carrier protein and preservative for stabilisation in the following formats:

<b>Concentrate:</b>	1 ml	(Cat. No. AIB-30018)
<b>Concentrate:</b>	0.5 ml	(Cat. No. AIB-30019)
<b>Ready-to-use:</b>	6 ml	(Cat. No. AIB-30017)

#### Dilution of primary antibody

Dilution of Nordic BioSites' concentrated antibody depends on the detection system used. The final working dilution must always be determined by the user. The elaboration of staining protocol should be done by an experienced specialist. For Nordic BioSites' recommendations see chapter 'Staining procedure'.

#### Storage and handling

The antibody should be stored at 2-8°C without further dilution.

Dilutions of the concentrated antibody should be done in a suitable antibody dilution buffer (e.g. BCB-20005/BCB-20006 from Nordic BioSite). The diluted antibody should be stored at 2-8°C after use. Stability of this working solution depends on various parameters and has to be confirmed by appropriate controls. The antibody provided is suitable for use until the expiry date indicated on the label, if stored at 2-8°C. Do not use product after the expiry date. Positive and negative controls should be run simultaneously with all specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Nordic BioSites' technical support or your local distributor.

#### Precautions

Use through qualified personnel only.

Wear protective clothing to avoid contact of reagents and specimens with eye, skin and mucous membranes. If reagents or specimens come in contact with sensitive area, wash with large amounts of water.

Microbial contamination of the reagent must be avoided, since otherwise non-specific staining may occur.

Sodium azide (NaN<sub>3</sub>), used for stabilisation, is not considered hazardous material in the concentration used. Reaction of sodium azide with lead or copper in drainage pipes can result in the formation of highly explosive metallic

azides. Sodium azide should be discarded in a large volume of running water to avoid formation of deposits. Material safety data sheets (MSDS) are available upon request.

### Staining procedure

Refer to the following table for conditions specifically recommended for this antibody. Also refer to detection system data sheets for guidance on specific staining protocols or other requirements.

Parameters	Nordic BioSite recommendations
*Pre-treatment	Heat Induced Epitope Retrieval in Citrate Buffer pH 6.0 (BCB-20015 / BCB-20016)
*Control tissue	Plasmacytoma or normal bone marrow
*Working dilution	1:200 (for concentrates)
*Incubation time	60 minutes

### Quality control

The recommended positive control tissues for this antibody are normal bone marrow or a plasmacytoma. We recommend carrying out a positive and a negative control with every staining run. Please refer to the instructions of the detection system for guidance on general quality control procedures.

### Troubleshooting

If you observe unusual staining or other deviations from the expected results please read these instructions carefully, refer to the instructions of the detection system for relevant information.

### Expected results

This antibody stains positive in the cytoplasmic membrane of plasma cells in formalin-fixed, paraffin-embedded tissue sections. Interpretation of the staining results is solely the responsibility of the user. Any experimental result should be confirmed by a medically established diagnostic procedure.

### Limitations of the Procedure

Immunohistochemistry is a complex technique involving both histological and immunological detection methods. Tissue processing and handling prior to immunostaining, for example variations in fixation and embedding or the inherent nature of the tissue can cause inconsistent results (Nadji and Morales, 1983). Endogenous peroxidase, alkaline phosphatase or biotin may cause non-specific staining depending on the detection system used. Tissues containing Hepatitis B Surface Antigen (HBsAg) may give false positive results with HRP (horse radish peroxidase) detection systems (Omata *et al*, 1980). Inadequate counterstaining and mounting can influence the interpretation of the results. Nordic BioSite warrants that the product will meet all requirements described from its shipping date until the expiry date is reached, if the product is stored and utilised as recommended. No additional guarantees can be given. Under no circumstances shall Nordic BioSite be liable for any damages arising out of the use of the reagent provided.

### Performance characteristics

Nordic BioSite has conducted studies to evaluate the performance of the antibody for use with a standard detection system. The product has been found to be sensitive and specific to the antigen of interest with minimal or no cross-reactivity.

### Bibliography

Sanderson RD *et al*. Cell Regul 1:27-35, 1989  
 Chilosi M *et al*. Mod Pathol 12:1101-1106, 1999  
 Bayer-Garner IB *et al*. Mod Pathol 14:1052-1058, 1999  
 Wei A and Juneja S. J Clin Pathol 56:406-411, 2003








Coupland SE *et al*. Br J Ophthalmol 89:352-359, 2005  
 Nadji M and Morales AR Ann N.Y. Acad Sci 420:134-9, 1983  
 Omata M *et al*. Am J Clin Pathol 73(5): 626-32, 1980

June 24, 2008

Rev: A0608

Doc: AIB-30017/30018/30019

Explanations of the symbols on the product label:

	Catalog Number Bestellnummer Reference du catalogue		Batch Code Chargenbezeichnung Code du lot	 Manufacturer  Nordic BioSite AB Propellervägen 4A S-183 62 Täby Sweden Tel: +46 (0)8 5444 33 40 Fax: +46 (0)8 756 94 90 info@nordicbiosite.com www.nordicbiosite.com
	Use By Verwendbar bis Utiliser jusqu'à		In Vitro Diagnostic Medical Device In vitro Diagnostikum Dispositif médical de diagnostic in vitro	
	Consult Instructions for use Gebrauchsanweisung beachten Consulter les instructions d'utilisation		Temperature Limitation Lagerungstemperatur Limites de température	