


Disposable Microtome Blades

- › S35
- › R35
- › N35
- › N35HR
- › A35
- › A22
- › S22
- › C35
- › S35L
- › S35LL
- › High Profile



pfm medical ag offers a broad range of FEATHER® disposable microtome blades for paraffin and cryostat sectioning for the whole range of individual applications in routine practice, research and industry.

www.pfmmedical.com



pfm medical ag offers a broad range of **FEATHER® disposable microtome blades** for paraffin and cryostat sectioning for the whole range of individual applications in routine practice, research and industry.

S35



Use

The S35 disposable microtome blade is a standard blade for routine practice and biopsies.

Details

- ▶ With excellent sectioning quality from the beginning
- ▶ Very long service life
- ▶ Short trimming time

Technical Data


- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

Ordering Information

Type S35

REF 207500000

R35



Use

The R35 disposable microtome blade is extra suitable for routine ribbon sectioning and for hard tissues.

Details

- ▶ Specially hardened cutting edge
- ▶ Short trimming time

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

Ordering Information

Type R35

REF 207500005

N35



Use

The N35 disposable microtome blade is particularly suitable for hard tissues.

Details

- ▶ Longer service life
- ▶ Suitable for all routine sectioning, also for soft tissue and biopsies
- ▶ Also suitable for ribbon sectioning after trimming
- ▶ Low compression

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

Ordering Information

Type N35

REF 207500006

N35HR



Use

The N35HR disposable microtome blade is particularly suitable for sectioning hard and difficult tissues.

Details

- ▶ Suitable for routine and ribbon sectioning
- ▶ Low compression
- ▶ Longer service life
- ▶ Short trimming time
- ▶ New coating
- ▶ Specially finished cutting edge

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

Ordering Information

Type N35HR

REF 207500014

A35



Use

The A35 disposable microtome blade allows the production of very thin sections and is suitable for hard and soft tissues.

Details

- ▶ Very long service life (also for hard tissue)
- ▶ Short trimming time with corresponding clearance angle
- ▶ Very sharp

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

Ordering Information

Type A35

REF 207500011

A22



Use

The A22 disposable microtome blade allows the production of very thin sections in routine use and research and is suitable for fibrous tissues (biopsies/fabrics).

Details

- ▶ Long service life
- ▶ Short trimming time with corresponding clearance angle
- ▶ Very sharp
- ▶ Also suitable for ribbon sectioning after trimming

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 22°
- ▶ **Dispenser:** 50 pcs.

Ordering Information

Type A22

REF 207500002

S22



Use

The S22 disposable microtome blade allows the production of extra-thin sections and is suitable for soft tissues/fabrics.

Details

- ▶ Long service life
- ▶ Short trimming time with corresponding clearance angle

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 22°
- ▶ **Dispenser:** 50 pcs.

Ordering Information

Type S22

REF 207500001

C35



Use

The C35 disposable microtome blade is specifically designed for cryostat sectioning.

Details

- ▶ Carbon steel blade
- ▶ Short trimming time

Technical Data

- ▶ **Material:** carbon steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

Ordering Information

Type C35

REF 207500003

S35L



Use

The S35L disposable microtome blade with a length of 120 mm is particularly suitable for large tissue blocks and routine sectioning.

Details

- ▶ Short trimming time

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 120 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

Ordering Information

Type S35L

REF 207500004

S35LL



Use

The extra-long S35LL disposable microtome blade with a length of 180 mm is suitable for large tissue blocks and routine sectioning.

Details

- ▶ Short trimming time

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 180 mm
- ▶ **Height:** 8 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 20 pcs.

Ordering Information

Type S35LL

REF 207500009

High Profile



Use

The High Profile disposable microtome blade is suitable for all paraffin embedded tissue blocks in routine applications.

Details

- ▶ Short trimming time

Technical Data

- ▶ **Material:** stainless steel
- ▶ **Length:** 80 mm
- ▶ **Height:** 14 mm
- ▶ **Blade angle:** 35°
- ▶ **Dispenser:** 50 pcs.

Ordering Information

Type High Profile

REF 207500010

Blade Depot



Use

The blade depot provides suitable storage for FEATHER® disposable microtome blades.

Details

- ▶ Storage capacity: approx. 200 blades
- ▶ Accessorie: protective cover blade deposit

Manufacturer

pfm medical ag
Wankelstraße 60
50996 Köln, Germany

Ordering Information

Blade Depot

REF 030000

Protective cover blade deposit

REF 030100

Ordering Information

Disposable Microtome Blades

Type	Colour coding	REF	Material	Length	Height	Blade angle	PU
S35	 dark blue	207500000	stainless steel	80 mm	8 mm	35°	50 pcs.
R35	 black	207500005	stainless steel	80 mm	8 mm	35°	50 pcs.
N35	 light green	207500006	stainless steel	80 mm	8 mm	35°	50 pcs.
N35HR	 dark green	207500014	stainless steel	80 mm	8 mm	35°	50 pcs.
A35	 light blue	207500011	stainless steel	80 mm	8 mm	35°	50 pcs.
A22	 pink	207500002	stainless steel	80 mm	8 mm	22°	50 pcs.
S22	 grey	207500001	stainless steel	80 mm	8 mm	22°	50 pcs.
C35	 white	207500003	carbon steel	80 mm	8 mm	35°	20 pcs.
S35L	 blue	207500004	stainless steel	120 mm	8 mm	35°	20 pcs.
S35LL	 green	207500009	stainless steel	180 mm	8 mm	35°	20 pcs.
High Profile	 brown	207500010	stainless steel	80 mm	14 mm	35°	50 pcs.

Blade Depot

REF	Capacity	PU
030000	200 blades	1 pc.

Protective cover blade deposit

REF	PU
030100	1 pc.

Contact

Should you have any questions our Customer Solutions Team is happy to advise you.

service@pfmmedical.com
T +49 (0)2236 9641-220
F +49 (0)2236 9641-51

www.pfmmedical.com - Your information source for the products of pfm medical ag.

pfmmedical

pfm medical ag
Wankelstraße 60
50996 Köln, Germany

Exclusive Distributor
for Europe

Certified by
DIN EN ISO 13485

Manufacturer
FEATHER® Safety Razor Co. Ltd.



Registered Trademark in Japan
Products are subject to alterations